



PHILIPS

AVENT

Improved flow rates for Philips Avent  
Natural 2.0 teats


innovation ✦ you

# Philips Avent Natural 2.0 supports even better the needs of growing babies thanks to the improved teat


- Ultra soft textured material (teats 0m+ -> 3m+)
- Spiral design for extra flexibility
- Improved flow rates
- Advanced anti colic valve



Natural 2.0 teats allow a natural tongue movement similar to breastfeeding, for a successful feeding



# The Natural 2.0 optimized flows are closer to breast flow rates



**Breastfeeding flow rates are considered “slow”, so slower flow rates facilitate the combination between breast and bottle**

When flow rates are too high:

- higher risk of teat rejection
- Possible breast rejection (“lazy baby”)
- distress signals, like coughing, gagging, spitting, spilling, etc.

# New Natural 2.0

0m+ 1m+ 3m+

Look and feel	Soft, silky texture spiral teat		
Ultrasoft silky flexible spiral teat			
Holes in teat	1	1	2
# of valves	1	1	1

6m+

Smooth	
Flexible bite resistant teat	
	3
	1

3m+

Smooth	
	--
	2

# Old Natural 1.0

0m+ 1m+ 3m+

Look and feel
Wide breast shaped teat
Holes in teat
# of valves

Smooth			
1	2	3	4
2	2	2	2

6m+

Smooth	
	--
	2

3m+

Printed copies are unaltered.





